














































## SYMBOLS GLOSSARY














Explanation of symbols used on the packaging of medical devices of ZARYS International Group sp. z o.o. sp.k.












	CE-mark
	CE-Mark with Notified Body number
	Manufacturer
	Authorized representative in the European Community
	Date of manufacture
	Use-by-date
	Batch code
	Reference number
	Importer
	Distributor
	Country of manufacture












	<p>Sterilized using ethylene oxide</p>
	<p>Sterilized using irradiation</p>
	<p>Sterilized using steam or dry heat</p>
	<p>Do not re-sterilize</p>
	<p>Non-sterile</p>
	<p>Do not use if package is damaged</p>
	<p>Single sterile barrier system</p>
	<p>Single sterile barrier system with protective packaging inside</p>
	<p>Single sterile barrier system with protective packaging outside</p>
	<p>Double sterile barrier system</p>
	<p>Fragile, handle with care</p>
	<p>Keep away from sunlight</p>












	<p>Keep dry</p>
	<p>Temperature limits to which the medical device can be safely exposed</p>
	<p>Humidity limitation</p>
	<p>Do not re-use</p>
	<p>Single patient multiple use</p>
	<p>Consult instructions for use</p>
	<p>Caution</p>
	<p>Contains or presence of natural rubber latex</p>
	<p>Not made with natural rubber latex</p>
	<p>Contains hazardous substances</p>

	<p>Contains or presence of 2-ethylhexyl phthalate (DEHP)</p>
	<p>Contains or presence of benzyl butyl phthalate (BBP)</p>
	<p>Contains or presence of dibutyl phthalate (DBP)</p>
	<p>Contains or presence of DEHP, BBP and DBP</p>
	<p>Phthalate-free</p>
	<p>Non-pyrogenic</p>
	<p>Drops per millilitre</p>
	<p>Liquid filter with pore size</p>
	<p>Gravity device</p>
	<p>One-way valve</p>
	<p>Flat angle</p>
	<p>Size</p>














	Medical Device
	Unique device identifier
	Packaging unit
	Packaging unit
	Inner diameter
	Outer diameter
	Length
	Length
	Cuff size
	Patient weight
	Maximum volume of air to inflate the cuff
	MRI safe
	Colour

	Size
	Grammage
	Non-woven type
	Number of threads in the gauze
	Number of plies
	Place of opening the package
	Place of opening the package
	Wash by hand
	Washing, normal process, maximum 40 degrees Celsius
	Washing, normal process, maximum 95 degrees Celsius
	Do not use bleach











	<p>Do not dry mechanically</p>
	<p>Do not iron</p>
	<p>Do not dry-clean</p>
	<p>Number of pairs</p>
	<p>Do not apply directly on wounds</p>
	<p>Contains superabsorbent</p>
	<p>In vitro diagnostic medical device</p>
	<p>Contains sufficient for &lt;n&gt; tests</p>
 <p>acceptable quality level EN 455-1 EN 374-2</p>	<p>Acceptable quality level EB 455-1, EN 374-2</p>
	<p>Powder-free gloves</p>
	<p>Number of gloves</p>









	<p>Number of pairs of gloves</p>
	<p>Contains or presence of Type IV allergen</p>
	<p>Stacking limit</p>
	<p>This side up</p>
	<p>Do not use hooks</p>
	<p>Do not use hooks</p>
	<p>Do not use sharp tools</p>
	<p>Do not trample</p>
	<p>Recyclable</p>
	<p>Recycling symbol: Low density polyethylene</p>
	<p>Recycling symbol: High density polyethylene</p>


















 <p>5 PP</p>	<p>Recycling symbol: Polypropylene</p>
 <p>20 PAP</p>	<p>Recycling symbol: Corrugated cardboard</p>
 <p>21 PAP</p>	<p>Recycling symbol: Flat cardboard</p>
 <p>22 PAP</p>	<p>Recycling symbol: Paper</p>
 <p>84 PAP+PE+ALU</p>	<p>Recycling symbol: Paper and cardboard/plastic/aluminium</p>
 <p>07 o</p>	<p>Recycling symbol: Other plastics</p>
	<p>Recycling symbol: Mixed</p>
	<p>Mark of conformity with the requirements of Ukrainian market</p>
 <p>UA.TR.116</p>	<p>Mark of conformity with the requirements of the Ukrainian market with the number of the conformity assessment body</p>
	<p>Mark of conformity with the requirements of UK market</p>
	<p>Certification according to ISO 14001</p>
	<p>Certification according to ISO 13485</p>
	<p>Authorised representative for the Swiss market</p>













Explanation of symbols used on the packaging of personal protective equipment of ZARYS International Group sp. z o.o. sp.k.










	CE-Mark
	CE-Mark with Notified Body number
 <p>Producent/ Manufacturer</p>	Manufacturer
 <p>Upoważniony przedstawiciel w Unii/ Authorized representative in the EU</p>	Authorized representative in the EU
 <p>Importer/ Importer</p>	Importer
 <p>Dystrybutor/ Distributor</p>	Distributor
 <p>Data produkcji/ Date of manufacture</p>	Date of manufacture
 <p>Użyć do daty/ Use-by-date</p>	Use-by-date
 <p>Kod partii/Batch code</p>	Batch code
 <p>Numer katalogowy/ Catalogue number</p>	Catalogue number










 <p>Trzymać z dala od światła słonecznego/ Keep away from sunlight</p>	Keep away from sunlight
 <p>Chronić przed wilgocią/ Keep dry</p>	Keep dry
 <p>Dopuszczalna temperatura/ Temperature limit</p>	Indicates the temperature limits to which the PPE can be safely exposed
 <p>Ograniczenie wilgotności/ Humidity limitation</p>	Humidity limitation
 <p>Nie używać powtórnie/ Do not re-use</p>	Do not re-use
 <p>Zawiera lateks/ Contains latex</p>	Contains or presence of natural rubber latex
 <p>Nie zawiera lateksu/ Latex-free</p>	Not made with natural rubber latex
 <p>Zajrzyj do instrukcji używania/ Consult instructions for use</p>	Consult instructions for use

 Rozmiar/ Size	Size
	Protection against chemicals
<b>TYP 4 B</b> <b>TYPE 4B</b>  <b>EN14605:2005</b> <b>+A1:2009</b>	Protection against spraying of liquid chemicals, type 4
<b>TYP 5 B</b> <b>TYPE 5B</b>  <b>EN ISO13982-1:2004</b> <b>+A1:2010</b>	Protection against airborne particles, type 5
<b>TYP 6 B</b> <b>TYPE 6B</b>  <b>EN13034:2005</b> <b>+ A1:2009</b>	Protection against liquid chemicals, type 6
 <b>EN1149-5:2008</b>	Protection against static electricity
 <b>EN14126:2003</b> <b>+ AC:2004</b>	Protection against micro-organisms

 <b>EN1073-2:2002</b>	<p>Protection against particulate radioactive contamination</p>
 Koniec okresu przechowywania/ End of shelf life	<p>End of shelf life</p>
 Patrz informacje dostarczone przez producenta/ See information supplied by the manufacturer	<p>See information supplied by the manufacturer</p>
 Zakres temperatur w warunkach przechowywania/ Temperature range of storage conditions	<p>Temperature range of storage conditions</p>
 Maksymalna wilgotność względna w warunkach przechowywania/ Maximum relative humidity of storage conditions	<p>Maximum relative humidity of storage conditions</p>
<p><b>EN ISO 374-1:2016 +A1:2018 TYPE X</b></p>  [(UVW)XYZ]	<p>Protection against chemicals</p>
<p><b>EN ISO 374-5: 2016</b></p>  <b>VIRUS</b>	<p>Protection against exposure to infectious agents</p>
 protection against blood-borne pathogens	<p>Protection against blood-borne pathogens</p>

 <p>protection against chemotherapy drugs</p>	<p>Protection against chemotherapy drugs</p>
 <p>acceptable quality level EN 455-1 EN 374-2</p>	<p>Acceptable quality level EN 455-1, EN 374-2</p>
	<p>Minimum glove length 240mm</p>
	<p>Acceptable quality level EN 455-1, EN 374-2</p>
	<p>Protection against pathogens</p>
	<p>Glove stretches in all directions</p>
	<p>Can be used with touch screens</p>
	<p>Contact with food</p>
	<p>Number of gloves</p>
	<p>Number of pairs of gloves</p>
 <p>Fit For Special Purpose (Palm protection)</p>	<p>See information supplied by the manufacturer</p>
	<p>See user information</p>

 <p>Nie prać/ Do not wash</p>	Do not wash
 <p>Nie używać wybielacza/ Do not bleach</p>	Do not bleach
 <p>Nie suszyć mechanicznie/ Do not tumble dry</p>	Do not tumble dry
 <p>Nie prasować/ Do not iron</p>	Do not iron
 <p>Nie czyścić chemicznie/ Do not dry clean</p>	Do not dry clean
 <p>Materiał łatwopalny/ Flammable fabric</p>	Flammable fabric
	Packaging unit
	Stacking limit
	This side up

	Do not use hooks
	Do not use sharp tools
	Do not trample
	Recyclable
	Recycling symbol: Low density polyethylene
	Recycling symbol: Corrugated cardboard
	Recycling symbol: Flat cardboard
	Recycling symbol: Paper
	Recycling symbol: Mixed